
— YOUR DNA KIT

Three steps. Seven minutes.

Welcome to MySpawn. Everything you need to collect and return your sample is in this box. Your legacy starts here.

WHAT'S IN THIS KIT



125 ml Sample Vial

Sterile, sealed, single-use



Return Envelope

Pre-addressed to the vault

INSTRUCTIONS

01 — 03

1

Collect your sample

Wash and dry your hands. Pull **5–10 strands of hair with the roots attached** (pulled, not cut — the root holds the DNA), and clip **3–4 fingernail trimmings**.

2

Seal them in the vial

Drop the hair and nail clippings straight into the vial. Twist the cap until it clicks. Keep it dry — no water or any liquid.

3

Mail it to the vault

Place the vial in the pre-addressed envelope — it routes to GenVault, our independent, FDA-registered, ISO 9001/20387 & CAP-accredited biorepository. Affix postage and post it.